

DB Antigen retrieval solution pH9 – 10x

For in vitro diagnostic use (IVD)

Cat. #: DB RETRIS-250 (250 ml) 10x concentrated
DB RETRIS-500 (500 ml) 10x concentrated

Product information:

Storage: 4-8 °C
Expiration: see on label
Preservative: 0.5% NaN₃
Buffer: 100 mM Tris-HCl, pH 9.0, 0.5% Tween-20, 10 mM EDTA

DB Antigen retrieval solution pH9 is designed for the application in immunohistochemistry to optimize the antibody performance with heat induced antigen retrieval on formalin-fixed, paraffin-embedded tissues.

Procedure:

Prepare a working-ready to use solution by diluting the concentrate DB RETRIS 1:10 in distilled water. The ready to use solution can be used for up to 3 months if stored at 4-8 °C.

Precautions:

1. For professional use only.
2. Wear appropriate protective equipment to avoid contact with eyes and skin.