

DB Antigen retrieval solution pH6 – 10x

For in vitro diagnostic use (IVD)

Cat. #: DB RETCIT-250 (250 ml) 10x concentrated
DB RETCIT-500 (500 ml) 10x concentrated

Product information:

Storage: 4-8 °C
Expiration: see on label
Preservative: 0.5% NaN₃
Buffer: 100 mM citrate buffer, pH 6, 0.5% Tween-20

DB Antigen retrieval solution pH6 is designed for the application in immunohistochemistry to optimize the antibody performance with heat induced antigen retrieval on formalin-fixed, paraffin-embedded tissues.

Procedure:

Prepare a working-ready to use solution by diluting the concentrated DB RETCIT 1:10 in distilled water. The ready to use solution can be used for up to 3 months if stored at 4-8 °C.

Precautions:

1. For professional use only.
2. Wear appropriate protective equipment to avoid contact with eyes and skin.