Anti - Erk 1,2

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KAT. ČÍSLO

DB 012-0.05 (50 µl)
DB 012-0.1 (100 µl)

WESTERN BLOT (WB) PROTOCOL

Roztoky pre Western blotting:

- Premývací pufo: 1x Tris Buffered Saline (TBS); 0,2% Triton X-100
- Blokovací pufo: 1x TBS; 0,2% Triton X-100; 8% sušené mlieko

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
IMUNOPRECIPITÁCIA (IP) PROTOKOL - NÁVOD NA POUŽITIE

1. Dilute sample (200 - 500 µg of total protein) with Correction Pufo (2.5% v/v Nonidet P-40, 5% w/v sodium deoxycholate, 0.5% w/v SDS in ddH₂O) in ratio 4:1 (v/v).
2. Add 5 µl of DB012 (anti-Erk 1,2, clone B20-U) rabbit clonal monospecific antibody, mix gently and incubate for 1 hour on ice.
3. Mix with 50 µl of ProteinG-Sepharose (washed with 10mM Tris-HCl, pH7.5*) and incubate for 30 minutes at 4°C with gentle shaking.
*NOTE: Washing of ProteinG-Sepharose with 10mM Tris-HCl, pH7.5: Resuspend 50 µl of ProteinG-Sepharose in 1 ml of 10mM Tris-HCl, pH7.5 by precise inverting of the tube several times. Centrifuge for 1 min, 900xg at 4°C and discard supernatant. Avoid of wasting ProteinG-Sepharose agarose gel beads during discarding. Repeat this procedure for 3 times.
4. Centrifuge ProteinG-Sepharose immunocomplex for 2 min, 900xg at 4°C and discard supernatant. Wash the pellet 3 times with 1 ml of RIPA Pufo (10mM TRIS-HCl, pH7.5, 140mM NaCl, 1% v/v Nonidet P-40, 1% w/v sodium deoxycholate, 0.5% w/v SDS in ddH₂O).
IMPORTANT: Avoid of wasting/discarding ProteinG-Sepharose immunocomplex.
5. Wash sediment with 1 ml of 10mM TRIS-HCl, pH7.5, centrifuge sample (agarose gel beads) for 1 min, 900xg at 4°C and discard the supernatant.
6. Dissociate immunocomplex from ProteinG-Sepharose with the help of Reduction Pufo (125mM TRIS-HCl, pH6.8, 3.3% SDS, 5% β-mercaptoethanol). Mix the sample with 30 µl of Reduction Pufo, shake gently and incubate for 5 minutes at 65°C.
7. Centrifuge at 3000xg for 5 minutes and transfer the supernatant (immunoprecipitated proteins) to new tube.
8. Separate the immunoprecipitated protein by 1D SDS-PAGE.

Representative picture of Erk 1,2 immunoprecipitated from HEK 293 cells, visualized with clonal rabbit anti- Erk 1,2 monospecific antibody (DB012) in Western blot analysis. Primary antibody dilution - 1:5,000.

lane 1 – molecular weight marker; lane 2 – positive control, Western blot of HEK 293 (100 µg) with anti-Erk 1, 2 rabbit clonal primary antibody; lane 3 – negative control, Immunoprecipitácia without primary antibody; lane 4 – immunoprecipitated Erk 1,2 (~40 and ~45 kDa proteins) from HEK 293 cells (500 µg of crude protein extract) using DB012 antibody

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