

DB detection system rabbit / mouse dual - HRP

For in vitro diagnostic use (IVD)

Cat. #: DB DET SYS (100 ml - 500 tests)

Product information

Storage: 2-8 °C

Expiration: see on label

DB rabbit/mouse dual-HRP detection system is a polymer based secondary antibody: HRP conjugate, designed for application in immunohistochemistry. It can be used with prediluted and concentrated antibodies (produced in rabbit and/or mouse) from different manufacturers.

Content:

HRP Rabbit/Mouse 100 ml ready-to-use solution containing secondary antibody conjugated with HRP, recognizing specifically rabbit and/or mouse IgG. Presented in buffered solution containing stabilizing protein and preservative.

Application:

After incubating the tissue section with primary antibody and appropriate washing procedure, apply the DB detection system-rabbit/mouse dual, ready to use for 25-30 minutes.