

DB primary antibody diluent - ready to use

For in vitro diagnostic use (IVD)

Cat. #: DB D-125 (125 ml)
DB D-250 (250 ml)

Product information:

Storage: 4-8 °C
Expiration: see on label
Preservative: 0.05% NaN₃
Buffer: 20 mM Tris-HCl, pH 8.0
Stabilizer: 10 mg/ml BSA

DB Antibody diluent is designed for the application in immunohistochemistry. The reagent is in a ready-to-use format and is used to dilute the concentrated primary antibodies. It is designed in order to optimize the antibody performance and reduce the possible nonspecific background staining.

Procedure:

Based on suggested dilutions, add the concentrated antibody directly to the diluent and mix thoroughly.